

Piezo Mini

MINI DENTAL PIEZO SURGERY

Minisize and lightweight
dual-screen design for smart dual control



MINI DESIGN

COMPACT & EFFICIENT

Balancing minimally invasive techniques & power output

5"

MAIN SCREEN



Full-touch interface for intuitive operation

1.6"

ROUND SCREEN



Real-time data display and adjust

OPTIMIZED DIMENSIONS

perfectly designed for a streamlined workspace

SLIM HANDPIECE

Operates at full power for 10 mins with a handpiece temperature of only 33°C, ensuring safety and comfort



Bone Mode



Endo Mode



Perio Mode



Custom Mode



Clean Mode

PARAMETER

VOLTAGE 100~240V 50/60Hz 120VA

FUSES 2×F10A 250V

FREQUENCY 26~38Hz

RIDAL AMPLITUDE 20-200μm

LATERAL AMPLITUDE ≤30μm

WATER FLOW 60-100ml

www.dentalgerate.com

PLUWA®
PLUS POWER
SN: S02524JUA1183
SG60P